



CHROMABOND Diamino

CHROMABOND® Diamino

QuEChERS method and pre-mixes

Within a few years after its development by Anastassiades et al. the QuEChERS method has gained a leading position for determination of pesticide residues in food by GC-MS or LC-MS allowing rapid and cheap clean-up of strong matrix-contaminated samples.

Standard clean-up of food samples

10 g sample are homogenised with 10 ml acetonitrile. After adding the internal standard the sample is shaken with 4 g MgSO₄ and 1 g NaCl and afterwards centrifuged. 1 ml of the supernatant is spiked with 25 mg CHROMABOND® Diamino and 150 mg MgSO₄ and shaken again. After centrifugation the supernatant is injected into the GC/MS.

Special silica phase for determination of pesticides in food samples

Base material silica, pore size 60Å, particle size 45µm, specific surface 500m²/g, pH stability 2 to 8. Primary and Secondary Amine functions (PSA), 5% C removes polar compounds (e. g. organic acids, pigments, sugars) from matrices like fruit or vegetables
similar phases: Supelclean PSA, Bond Elut PSA.

Recommended application: special SPE phase for quick and cheap determination of pesticides in strongly matrix-contaminated samples by GC. QuEChERS method = Quick Easy Cheap Effective Rugged Safe.

QuEChERS clean-up mixes containing 0.15 g CHROMABOND® Diamino each.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

