OPERATION MANUAL

LLG-uniSTIRRER 3 pro

Thank you very much for purchasing this
LLG-uniSTIRRER 3 pro.

Please read the following instructions carefully
and always keep this manual within easy reach.

Contents

1. Preface/Service	3	8. Stirring	12
2. Safety Instructions	3	9. Faults	12
3. Proper Use	5	10. Maintenance and Cleaning	13
4. Inspection	6	11. Associated Standards and Regulations	14
4.1 Receiving Inspection	6	12. Specifications	15
4.2 Listing of Items	6	13. Ordering Information	16
5. Control	7	14. Warranty	17
5.1 Control elements	7		
5.2 Display	8		
6. Trial Run	9		
7. Function: Heating	10		
7.1 Working with external temperature sensor	10		
7.2 Residual heat warning (HOT warning)	11		

1. Preface/Service

Thank you very much for purchasing the LLG-uniSTIRRER 3 pro. You should read this Manual carefully, follow the instructions and procedures, and be aware of all the cautions when using this instrument.

When help needed, you can always contact your supplier for technical support.

2. Safety Instructions

	<p>Warning!</p> <ul style="list-style-type: none"> • Read the operating instructions carefully before use. • Ensure that only trained staff works with the instrument.
	<p>Risk of burn!</p> <ul style="list-style-type: none"> • Caution when touch the housing parts and the hotplate (heating model) which can reach temperature of 280 °C. • Pay attention to the residual heat after switching off.
	<p>Protective ground contact!</p> <ul style="list-style-type: none"> • Make sure that socket must be grounded (protective ground contact) before use.

- When working with the stirrer please wear personal safety protection to avoid the risk of:
 - Splashing and evaporation of liquids
 - Release of toxic or combustible gases
- Set up the instrument in a spacious area on a stable, clean, non-slip, dry and fireproof surface. Do not operate the instrument in explosive atmospheres, with hazardous substances or under water.
- Gradually increase the speed, reduce the speed if:
 - Stirring bar breaks away due to high speed
 - The instrument is not running smoothly or container moves on the base plate
- Temperature must always be set at least 50 °C lower than the fire point of the media used.
- Be aware of hazards due to:
 - Flammable materials or media with a low boiling temperature
 - Overfilling of media
 - Unsafe container
- Process pathogenic materials only in closed vessels.
- Check the instrument and accessories prior to each use. Do not use damaged components. Safe operation is only guaranteed with the accessories described in the "Accessories" chapter. Accessories must be securely attached to the device and cannot come off by themselves. Always disconnect the plug before fitting accessories.
- When the external temperature sensor is needed, the tip of the measuring sensor must be at least 5-10 mm from vessel bottom and wall.
- The instrument can only be disconnected from the main power supply by pulling out the main or the connector plug.
- The voltage stated on the label must correspond to the main power supply.
- Ensure that the main power supply cable does not touch the hotplate. Do not cover the device.
- The instrument may only be opened by experts.
- Keep away from high magnetic field.

3. Proper Use

The LLG-uniSTIRRER 3 pro is designed for mixing and/or heating liquids in schools, laboratories or factories.

- Observe the minimum distances between the devices, between the device and the wall and above the assembly (min. 100 mm)


Figure 1

This device is not suitable for using in residential areas or other constraints mentioned in Chapter 1.

4. Inspection

4.1 Receiving Inspection

Unpack the equipment carefully and check for any damages, which may have arisen during transport. Please contact your supplier for technical support.

Note:

If there is any apparent damage to the system, please do not plug it into the power line.

4.2 Listing of Items

The package includes the following items:

Items	Qty
LLG-uniSTIRRER 3 pro	1
Power cable	1
User Manual	1
Temperature sensor PT1000	1
Support Rod	1
Clamp for PT1000	1

Table 1


Figure 2

5. Control

5.1 Control elements


Items	Descriptions
Temperature Control Knob Heat	Set the temperature parameters. The function "heating" is switched ON or OFF via push ON/OFF knob.
Speed Control Knob	Set the rotary speed. The function "Stirring" is switched ON or OFF via push ON/OFF knob.
LED Display	LED displays the real working state and all settings.
Probe	When the external temperature sensor PT1000 is plugged in, probe icon  is lit.
Power Switch	Switch ON or OFF the instrument.

Table 2

Power switch

Figure 3

5.2 Display


Figure 4

Characters	Descriptions
Temperature display area	When heating function is switched ON, the LED displays the temperature setting value and shifts to real value in 5 seconds. When the heating function is switched OFF and the hotplate temperature is still above 50 °C, LED displays "Hot", otherwise LED displays OFF.
Speed display area	When stirring function is switched ON, LED displays the speed setting value and flashes. The setting value does not flash until the current speed reaches the set

Table 3

6. Trial Run

- Make sure the required operating voltage and powersupply voltage match.
- Ensure the socket must be properly grounded.
- Plug in the power cable, ensure the power is on and begin initializing.
- Add the medium into the vessel with an appropriate stirring bar.
- Place vessel on the work plate.
- Set the target stirring speed and begin.
- Observe the stirring bar and LED display (digital model).
- Set the target temperature and start heating.
- Observe the LED display (digital hotplate model).
- Stop the heating and stirring functions.

If these operations above are normal, the device is ready to operate. If not, the device may be damaged during transportation, please contact your supplier for technical support.

Warning!

Forbid to transfer the vessel when the instrument working.

7. Function: Heating

The device is controlled by digital temperature control technology, which has two separate safe circuits. The hotplate is kept at a constant temperature by a digital control circuit. The hotplate temperature can also be monitored from a separate, adjustable safe circuit. The two temperature sensors internal for temperature control are built into the hotplate. The single external PT1000 can monitor the temperature of sample.

- Plug in the external PT1000.
- Set the temperature via rotating the temperature control knob slowly to the target value.
- When the heating function is switched on, the LED displays the temperature value on the left-hand side.
- The heating function is switched on or off by pushing heating knob.

The LLG-uniSTIRRER 3 pro automatically displays the last running speed and temperature parameters once turned on. To ensure the accuracy of the temperature inside the container, please use the external temperature sensor PT1000

7.1 Working with external temperature sensor


Figure 5

The external temperature sensor PT1000 is the standard accessory. When PT1000 sensor is connected and the temperature control knob is used, LED displays the temperature setting value and shifts to real value in 5 seconds. Safe circuit controls hotplate temperature under 320 °C. Comparing with the temperature control of the hotplate, the external temperature sensor can control the medium's temperature more precise. The heating function will be stopped automatically under abnormal conditions. Please operate follow the instructions below:

- Switch OFF the instrument.
- Ensure the external temperature sensor is inserted in the media heated.
- Switch ON the instrument and run heating function.

If the heating function does not work, please contact your supplier for technical support.

7.2 Residual heat warning (HOT)

In order to prevent the risk of burns from a hotplate, the hotplate stirrer has a residual heat warning function. When the heating function is switched off and the heating plate temperature is still above 50 °C, "Hot" will flash to warn that there is a hazard of burns from the hotplate. When the unit is powered off, the LED screen displays the temperature of hotplate and Hot alternately. When the hotplate temperature drops to below 50 °C, the unit will automatically switch off. If you want to turn off the LED immediately, just pull out the plug directly. When the plug is pulled out, the residual heat warning function cannot be run.

8. Stirring

The function "stirring" of the LLG-uniSTIRRER 3 pro is switched on or off via pushing on/off speed control knob. The speed range is 200 to 1500 rpm, in steps of 10 rpm.

9. Faults

- Instruments cannot be powered ON
 - Check whether the power line is unplugged
 - Check whether the fuse is broken or loose
- Fault in power on self test
 - Switch OFF the unit, then switch ON and reset the instruments to factory default setting.
- Stir speed cannot reach set point
 - Excessive medium viscosity may cause abnormal speed reduction of the motor
- Unit cannot be powered off when switched off.
 - Check if the residual heat warning function is still ON and hotplate temperature is above 50 °C (the LED screen still work and "Hot" flash).

If these faults are not resolved, please contact your supplier.

10. Maintenance and Cleaning

- Proper maintenance can keep instruments working properly and lengthen its lifetime.
- Do not spray cleanser into the instrument when cleaning.
- Unplug the power line when cleaning.
- Use recommended cleansers only:

Dyes	Isopropyl alcohol
Construction materials	Water containing tenside / Isopropyl alcohol
Cosmetics	Water containing tenside / Isopropyl alcoho
Foodstuffs	Water containing tenside
Fuels	Water containing tenside

Table 4

- Wear the proper protective gloves during cleaning of the instrument.
- If you require maintenance service, you must clean the instrument in advance to avoid pollution of hazardous substances. Please send back the device in the original packing.
- If the instrument will not use for a long time, please switch off and place in a dry, clean, room temperature and stable location.

Note:

Electronic devices do not clean with cleanser.

11. Associated Standards and Regulations

Construction in accordance with the following safety standards:

EN 61010-1

UL 3101-1

CAN/CSA C22.2(1010-1) EN 61010-2-10

Construction in accordance with the following EMC standards:

EN 61326-1

Associated EU guidelines:

EMC-guidelines: 89/336

EWG Instrument guidelines: 73/023/EWG

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note:

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

12. Specifications

Items	Specifications LLG-uniSTIRRER 3 pro	Items	Specifications LLG-uniSTIRRER 3 pro
Residual heat warning	50 °C	Temperature display accuracy [°C]	1
Dimensions (mm)	220 × 160 × 95	Safety temperature [°C]	320
Weight [kg]	1.4	External Temperature sensor	PT1000
Permitted ambient temperature [°C]	5-40	Accuracy of heating temperature with PT1000 [°C]	± 0,5
Permitted relative humidity [%]	80 %	Voltage [VAC]	100~120/200~240
Protection class acc. to DIN 60529	*IP21	Frequency [Hz]	50/60
Rotary speed display	LED	Power [W]	515/15
Plate material	Ceramic coated stainless steel	Stirring positions	1
Dimensions of workplate (mm)	φ135	Max. stirring quantity (H ₂ O) [l]	3
Heating power [W]	500	Max. magnetic bar [L×Ø, mm]	50
Temperature range [°C]	RT-280, increment: 1	Motor type	DC motor
Temperature display	LED	Max. power input of motor [W]	5
		Max. power output of motor [W]	3
		Speed range [rpm]	200-1500, increment:10

Table 5

13. Ordering Information

Cat No.	Descriptions
6.263 400	LLG-uni <i>STIRRER</i> 3 pro with EU plug
6.263 401	LLG-uni <i>STIRRER</i> 3 pro with UK plug
6.263 415	Temperature sensor for LLG-uni <i>STIRRER</i> 3 pro, length 230 mm, sensor length 198 mm
6.263 417	Temperature sensor glass coated for LLG-uni <i>STIRRER</i> 3 pro, length 230 mm, sensor length 198 mm
6.263 416	Rod and support clamp of PT1000

14. Warranty

This product is warranted to be free from defects in material and workmanship for a period of three (3) years from date of purchase.

This warranty is valid only if the product is used for its intended purpose and within the guidelines specified in this instruction manual. This warranty does not cover damage caused by accident, neglect, misuse, improper service, natural forces or other causes not arising from defects in original material or workmanship. This warranty does not cover any incidental or consequential damages, commercial loss or any other damages from the use of this product.

The warranty does not cover damage to paint or finish and defects or damages caused by physical and chemical abuse or normal wear and tear. The warranty is invalidated by any nonfactory modification, which will immediately terminate all liabilities on us for the products or damages caused by its use. The buyer and its customer shall be responsible for the product or use of products as well as any supervision required for safety. If requested the products must be returned to the distributor in well packed and insured manner and all shipping charges must be paid.

Information on the disposal of electrical and electronic devices in the European Community:

The disposal of electrical devices is regulated within the European Community by national regulations based on the European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). According to these regulations, any devices supplied after 13.06.2005 in the business to business sphere, to which this product is assigned, may no longer be disposed of in municipal or domestic waste. They are marked with the following symbol to indicate this.

As disposal regulations within the EU may vary from country to country, please contact your supplier if necessary.

PRODUCT DISPOSAL

In case the product is to be disposed of, the relevant legal regulations are to be followed.

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Lab Logistics Group GmbH