



Jeringas para el automuestreador CG de Shimadzu

Trajan Scientific Europe Ltd

Jeringas para el automuestreador CG de Shimadzu

Para Shimadzu AOC-14, AOC-17, AOC-20 y AOC-20i. Punta de la aguja: cónica.



Code	Description	Packaging
LLG06205990	syringe 5F-S-0.63 5 µl	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

