

**Microtubos Safe-Lock, PP****Eppendorf SE****Microtubos Safe-Lock, PP**

- Un pequeño gancho del cierre del tapón se engancha en el reborde del microtubo
- Autoclavables estando abiertos (121 °C, 20 min.)
- Se pueden abrir y cerrar fácilmente con una mano
- Con superficie esmerilada de rotulación
- El gancho del tapón impide que salte por accidente (por ej. durante la termostatización)
- Máxima estabilidad mecánica para la centrifugación
- Graduados
- Los recipientes de reacción Safe-lock también están disponibles en calidad Eppendorf Biopur®
- Los recipientes de reacción con cierre de seguridad de 0,5 ml y 1,5 ml de capacidad pueden ser centrifugados hasta 30 000 x g



| Code        | Description                              | Packaging |
|-------------|------------------------------------------|-----------|
| LLG09409320 | SAFE-LOCK reaction vessels0.5 ml,co/500  | 500 pz.   |
| LLG09409331 | SAFE-LOCK reaction vessels1.5 ml,c/ 1000 | 1000 pz.  |
| LLG09409332 | SAFE-LOCK reaction vessels,cap. 1./ 1000 | 1000 pz.  |
| LLG09409333 | SAFE-LOCK reaction vessels,cap. 1./ 1000 | 1000 pz.  |
| LLG09409334 | SAFE-LOCK reaction vessels,cap. 1./ 1000 | 1000 pz.  |
| LLG09409335 | SAFE-LOCK reaction vessels,cap. 1./ 1000 | 1000 pz.  |
| LLG09409336 | SAFE-LOCK reaction vessels,cap. 1./ 1000 | 1000 pz.  |
| LLG09409341 | SAFE-LOCK reaction vessels2 ml,col/ 1000 | 1000 pz.  |
| LLG09409342 | SAFE-LOCK reaction vessels,cap. 2 / 1000 | 1000 pz.  |
| LLG09409343 | SAFE-LOCK reaction vessels,cap. 2 / 1000 | 1000 pz.  |
| LLG09409344 | SAFE-LOCK reaction vessels,cap. 2 / 1000 | 1000 pz.  |
| LLG09409345 | SAFE-LOCK reaction vessels,cap. 2 / 1000 | 1000 pz.  |
| LLG09409346 | SAFE-LOCK reaction vessels,cap. 2 / 1000 | 1000 pz.  |
| LLG09409352 | SAFE-LOCK reaction vessels,cap. 0.5/500  | 500 pz.   |
| LLG09409353 | SAFE-LOCK reaction vessels,cap. 0.5/500  | 500 pz.   |
| LLG09409354 | SAFE-LOCK reaction vessels,cap. 0.5/500  | 500 pz.   |
| LLG09409355 | SAFE-LOCK reaction vessels,cap. 0.5/500  | 500 pz.   |
| LLG09409356 | SAFE-LOCK reaction vessels,cap. 0.5/500  | 500 pz.   |
| LLG09409357 | Safe-Lock reaction vessel 0,5 mlass/500  | 500 pz.   |
| LLG09409358 | Safe-Lock reaction vessels 1,5 mla/ 1000 | 1000 pz.  |



| Code        | Description                              | Packaging |
|-------------|------------------------------------------|-----------|
| LLG09409359 | Safe-Lock reaction vessels 2,0 mla/ 1000 | 1000 pz.  |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6ª Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

