



Homogeneizador DUALL[®] con pistilo de vidrio

Homogeneizador DUALL[®] con pistilo de vidrio

Con área cónica y cilíndrica para triturar y homogeneizar uniformemente tejidos de forma efectiva. Especialmente indicado para tejidos celulares como músculos, corazón y pulmones. El pistilo es suficientemente estable para funcionar con motor. La distancia entre el pistilo y el mortero es de 0,1 a 0,15 mm. Se suministran repuestos, todos intercambiables.



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09651637 | DUALL homogenizers with glass pestle1ml | 1 pz |
| LLG09651638 | DUALL homogenizers, with glass pestle3m | 1 pz |
| LLG09651639 | Duall homogenizers,with glass pestle5 m | 1 pz |
| LLG09651641 | DUALL homogenizers with glass pestle15m | 1 pz |
| LLG09651642 | DUALL homogenizers with glass pestle30m | 1 pz |
| LLG09651643 | DUALL homogenizers with glass pestle50m | 1 pz |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

