



Homogeneizador DUALL[®] con pistilo de vidrio

Homogeneizador DUALL[®] con pistilo de vidrio

Con área cónica y cilíndrica para triturar y homogeneizar uniformemente tejidos de forma efectiva. Especialmente indicado para tejidos celulares como músculos, corazón y pulmones. El pistilo es suficientemente estable para funcionar con motor. La distancia entre el pistilo y el mortero es de 0,1 a 0,15 mm. Se suministran repuestos, todos intercambiables.



Code	Description	Packaging
LLG09651637	DUALL homogenizers with glass pestle1ml	1 pz
LLG09651638	DUALL homogenizers, with glass pestle3m	1 pz
LLG09651639	Duall homogenizers,with glass pestle5 m	1 pz
LLG09651641	DUALL homogenizers with glass pestle15m	1 pz
LLG09651642	DUALL homogenizers with glass pestle30m	1 pz
LLG09651643	DUALL homogenizers with glass pestle50m	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

