

**Separating funnels acc. to Squibb, borosilicate glass 3.3****Lenz****Separating funnels acc. to Squibb, borosilicate glass 3.3**

With PTFE plug. Available with or without scale to read the approximate volume.

With PP-stopper .



Code	Description	Packaging
LLG09203322	Separating funnels acc. to Squibb, borosilicate glass 3.3, without scale, 50 ml, Bore 2.5 mm, Socket 19/26 NS	1 pz.
LLG09203323	Separating funnels acc. to Squibb, borosilicate glass 3.3, without scale, 100 ml, Bore 2.5 mm, Socket 19/26 NS	1 pz.
LLG09203372	Separating funnels acc. to Squibb, borosilicate glass 3.3, with scale, 50 ml, Bore 2.5 mm, Socket 19/26 NS	1 pz.
LLG09203373	Separating funnels acc. to Squibb, borosilicate glass 3.3, with scale, 100 ml, Bore 2.5 mm, Socket 19/26 NS	1 pz.
LLG09203325	Separating funnels acc. to Squibb, borosilicate glass 3.3, without scale, 250 ml, Bore 4 mm, Socket 29/32 NS	1 pz.
LLG09203328	Separating funnels acc. to Squibb, borosilicate glass 3.3, without scale, 500 ml, Bore 4 mm, Socket 29/32 NS	1 pz.
LLG09203375	Separating funnels acc. to Squibb, borosilicate glass 3.3, with scale, 250 ml, Bore 4 mm, Socket 29/32 NS	1 pz.
LLG09203376	Separating funnels acc. to Squibb, borosilicate glass 3.3, with scale, 500 ml, Bore 4 mm, Socket 29/32 NS	1 pz.
LLG09203330	Separating funnels acc. to Squibb, borosilicate glass 3.3, without scale, 1000 ml, Bore 6 mm, Socket 29/32 NS	1 pz.
LLG09203332	Separating funnels acc. to Squibb, borosilicate glass 3.3, without scale, 2000 ml, Bore 6 mm, Socket 29/32 NS	1 pz.
LLG09203377	Separating funnels acc. to Squibb, borosilicate glass 3.3, with scale, 1000 ml, Bore 6 mm, Socket 29/32 NS	1 pz.
LLG09203378	Separating funnels acc. to Squibb, borosilicate glass 3.3, with scale, 2000 ml, Bore 6 mm, Socket 29/32 NS	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

