

**Laboratory threaded joints - GL 14 / 18 / 25****BOLA****Laboratory threaded joints - GL 14 / 18 / 25**

Enables GL threads to clamp easily to fluoroplastic or thick-walled tubing made out of other plastics, glass or metal. (Ideally suited for connection with Bola GL fittings, Bola pressure valves, Bola GL taps and Bola GL ball valves).

The laboratory threaded joints consist of a red screw cap made of glass-fibre reinforced ETFE, a sealing disc and a tapered sealing ring made of PTFE, as well as a V-ring made of glass-fibre reinforced ETFE. This patented system withstands pressures up to approx. 10 bar, is suitable for a vacuum (at room temperature in the case of sizes 14 to GL25) and is easy to assemble without tools. The substance that passes through only comes in contact with PTFE.

Bola laboratory threaded joints can, for example, be used to hold probes (including stainless steel probes), sensors, laboratory taps, submerged tubes, thermometers or as a cable feed.



Code	Description	Packaging
LLG09209801	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 0.8 mm, Neck thread 14 GL	
LLG09209802	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 1.0 mm, Neck thread 14 GL	1 pz.
LLG09209803	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 1.6 mm, Neck thread 14 GL	1 pz.
LLG09209813	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 1.6 mm, Neck thread 18 GL	
LLG09209804	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 2.0 mm, Neck thread 14 GL	
LLG09209814	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 2.0 mm, Neck thread 18 GL	
LLG09209805	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 2.4 mm, Neck thread 14 GL	
LLG09209806	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 3.0 mm, Neck thread 14 GL	1 pz.
LLG09209816	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 3.0 mm, Neck thread 18 GL	1 pz.
LLG09209836	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 3.0 mm, Neck thread 25 GL	
LLG09209807	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 3.2 mm, Neck thread 14 GL	1 pz.
LLG09209817	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 3.2 mm, Neck thread 18 GL	1 pz.
LLG09209837	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 3.2 mm, Neck thread 25 GL	
LLG09209808	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 4.0 mm, Neck thread 14 GL	1 pz.
LLG09209818	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 4.0 mm, Neck thread 18 GL	1 pz.
LLG09209838	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 4.0 mm, Neck thread 25 GL	
LLG09209809	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 6.0 mm, Neck thread 14 GL	1 pz.
LLG09209819	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 6.0 mm, Neck thread 18 GL	1 pz.
LLG09209839	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 6.0 mm, Neck thread 25 GL	1 pz.
LLG09209810	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 6.35 mm, Neck thread 14 GL	1 pz.



Code	Description	Packaging
LLG09209820	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 8.0 mm, Neck thread 18 GL	1 pz.
LLG09209840	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 8.0 mm, Neck thread 25 GL	1 pz.
LLG09209821	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 10.0 mm, Neck thread 18 GL	1 pz.
LLG09209841	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 10.0 mm, Neck thread 25 GL	1 pz.
LLG09209842	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 12.0 mm, Neck thread 25 GL	1 pz.
LLG09209843	Laboratory threaded joints - GL 14 / 18 / 25, For external tube diam. 14.0 mm, Neck thread 25 GL	



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

