



Guard columns for EC columns NUCLEODUR Spinx RP, 5 μ m

Guard columns for EC columns NUCLEODUR® Spinx RP, 5 μ m

distinct selectivity based on bifunctional surface coverage - USP L1 and USP L11

pore size 110 Å, particle sizes 1,8 μ m, 3 μ m and 5 μ m; 14 %C

high density of covalently bonded silanes for tailing-free peaks widens the scope for method development

pH stability 1 - 10

suitable for **LC/MS** due low bleeding characteristics

high reproducibility and consistent quality due to tight QC procedures

range of application: quinolone antibiotics, sulfonamides, xanthines, substituted aromatics





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

