



## Guard columns for EC columns NUCLEODUR 100-5 CN, 5 $\mu$ m

### Guard columns for EC columns NUCLEODUR® 100-5 CN, 5 $\mu$ m

pore size 110 Å, particle sizes 3 $\mu$ m and 5 $\mu$ m; 7%C - USP L10

multi-mode columns (RP and NP)

widens the scope in selectivity

different retention characteristics compared to C8 and C18

stable against hydrolysis at low pH values, working range pH 1-8

high reproducibility from lot to lot

ideal for: tricyclic antidepressants, steroids, organic acids





**CARLO ERBA Reagents S.A.S.**  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

