



Tube tests NANOCOLOR Part 2

Tube tests **NANOCOLOR®** Part 2

Using **NANOCOLOR®** bar coded test tubes, the measurement is automatically performed after inserting the cuvette into the cuvette slot. The contact with any dangerous chemicals is avoided because the measurement is performed directly in the test tube. This guarantees maximum safety for the user.



Code	Description	Packaging
LLG09304622	Tube testsNANOCOLOR® Part 2, Type Nitrite 2	
LLG04005405	Tube testsNANOCOLOR® Part 2, Type Peroxide 2*	
LLG04005404	Tube testsNANOCOLOR® Part 2, Type Methanol 15	
LLG04005396	Tube testsNANOCOLOR® Part 2, Type Silver 3	
LLG04005403	Tube testsNANOCOLOR® Part 2, Type Thiocyanate 50	
LLG04005399	Tube testsNANOCOLOR® Part 2, Type Molybdenum 40	
LLG04005400	Tube testsNANOCOLOR® Part 2, Type POC 200	
LLG04005397	Tube testsNANOCOLOR® Part 2, Type Organic acids 3000	
LLG04005402	Tube testsNANOCOLOR® Part 2, Type Starch 100	
LLG04005406	Tube testsNANOCOLOR® Part 2, Type TTC / Sludge activity 150**	
LLG04005398	Tube testsNANOCOLOR® Part 2, Type ortho- and total Phosphate 45	
LLG06242072	Tube testsNANOCOLOR® Part 2, Type Anionic surfactants 4	20 test
LLG09304623	Tube testsNANOCOLOR® Part 2, Type ortho- and total Phosphate 1	19 pz
LLG06231412	Tube testsNANOCOLOR® Part 2, Type Residual hardness 1	
LLG09921213	Tube testsNANOCOLOR® Part 2, Type Zinc 4	





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

