

**EPR sample tubes, quartz glass****Hilgenberg****EPR sample tubes, quartz glass**

These EPR tubes are manufactured from upmarket quartz glass with only slightest amounts of paramagnetic trace elements. The tubes are delivered single packed and sealed in plastic foil.



| Code        | Description  | Packaging |
|-------------|--|-----------|
| LLG06285951 | EPR sample tubes, quartz glass, Outer diam. 3.95 mm, Int. diam. 3.35 mm, Wall thickness 0.3 mm, Length 250 mm  | 10 pz.    |
| LLG06285950 | EPR sample tubes, quartz glass, Outer diam. 2.95 mm, Int. diam. 2.35 mm, Wall thickness 0.3 mm, Length 250 mm  | 10 pz.    |
| LLG06285952 | EPR sample tubes, quartz glass, Outer diam. 4.95 mm, Int. diam. 4.19 mm, Wall thickness 0.38 mm, Length 250 mm | 10 pz.    |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

