


**EPR sample tubes, quartz glass****Hilgenberg****EPR sample tubes, quartz glass**

These EPR tubes are manufactured from upmarket quartz glass with only slightest amounts of paramagnetic trace elements. The tubes are delivered single packed and sealed in plastic foil.



Code	Description	Packaging
LLG06285950	EPR sample tubes, quartz glass, Outer diam. 2.95 mm, Int. diam. 2.35 mm, Wall thickness 0.3 mm, Length 250 mm	10 pz
LLG06285951	EPR sample tubes, quartz glass, Outer diam. 3.95 mm, Int. diam. 3.35 mm, Wall thickness 0.3 mm, Length 250 mm	10 pz
LLG06285952	EPR sample tubes, quartz glass, Outer diam. 4.95 mm, Int. diam. 4.19 mm, Wall thickness 0.38 mm, Length 250 mm	10 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

