

**Bottles (Karlsruhe bottles) with stoppers**

behr

**Bottles (Karlsruhe bottles) with stoppers**

The special bottle neck of the Karlsruhe bottles absorbs the water which is displaced by the stopper of the electrode. This construction guarantees clean working and ensures at the same time that the test results are not adulterated by unwanted entry of oxygen into the samples.



Code	Description	Packaging
LLG09920515	Bottles (Karlsruhe bottles) with stoppers, 100 ml, Length* 60	1 pz.
LLG09920514	Bottles (Karlsruhe bottles) with stoppers, 100 ml, Length* 20	1 pz.
LLG09920509	Bottles (Karlsruhe bottles) with stoppers, 250 ml, Length* 60	1 pz.
LLG09920513	Bottles (Karlsruhe bottles) with stoppers, 250 ml, Length* 20	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

