

**Accessories for shakers HS/KS 260****IKA****Accessories for shakers HS/KS 260****AS 260.1 Universal attachment**

Able to hold most types of container using adjustable tension rollers. Surface: 320 x 320mm

AS 260.2 Fixing clip platform

For shaking round bottom vessels, volumetric and Erlenmeyer flasks. Supplied without clips.

Max. number of clips (capacity held): 56 x AS 2.1 (25ml), 23 x AS 2.2 (50ml), 23 x AS 2.3 (100ml), 11 x AS 2.4 (250ml), 9 x AS 2.5 (500ml), 5 x AS 2.6 (1000ml)

AS 260.3 Petri dish attachment

For gentle mixing and suspension of cell cultures or growth media in petri dishes, culture bottles and containers. With integral anti-slide mat (PP).

AS 260.5 Separating funnel attachment

For HS 260 basic and control

For shaking, elution, extraction, gassing out, dissolving, enriching. The 3 tension rollers supplied are height-adjustable to accommodate different separating funnel sizes. The separating funnels are secured by means of O-rings (6 O-rings supplied).

Capacity: Number of funnels (pear shaped) x volume 6 x 50ml, 5 x 100ml, 3 x 250ml, 3 x 500ml



Code	Description	Packaging
LLG09838192	Accessories for shakers HS/KS 260, Type AS 260.3, Petri dish attachment, dim. (W x D x H) 410 x 370 x 32 mm	1 pz.
LLG09838191	Accessories for shakers HS/KS 260, Type AS 260.2, Gridded clip platform attachment*, dim. (W x D x H) 330 x 330 x 24 mm	1 pz.
LLG09838193	Accessories for shakers HS/KS 260, Type AS 260.5, Separating funnel attachment, dim. (W x D x H) 334 x 425 x 145 mm	1 pz.
LLG09838080	Accessories for shakers HS/KS 260, Type AS 260.1, Universal attachment, dim. (W x D x H) 425 x 335 x 135 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

