

**Selectron filter holders FP 025/1 and FP 050/1****Whatman****Selectron® filter holders FP 025/1 and FP 050/1**

In polysulphone for use with syringes for analytical and preparative use.
Filtration-media must be selected separately.

Applications of FP 025/1:

Clarification and sterile filtration of aqueous liquids and gases.

Applications of FP 050/1:

Sampling and flow rate applications. Ideal for filtrations using peristaltic pumps and for gas sterilisation.

Max. pressure:

Max. temperature: 7 bar

180 °C



Code	Description	Packaging
LLG09603181	Selectron® filter holders FP 025/1 and FP 050/1, Type FP 025/1, Diam. 25 mm	10 pz.
LLG09603185	Selectron® filter holders FP 025/1 and FP 050/1, Type FP 050/0, Diam. 50 mm	
LLG09603205	Selectron® filter holders FP 025/1 and FP 050/1, Type FP 050/1, Diam. 50 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

