



## Software EasySieve for Sieve Shakers AS control / AS 200 tap / AS jet

### Software EasySieve® for Sieve Shakers AS control / AS 200 tap / AS jet

Software for particle size analyses.

#### Features:

- Automatic registration, evaluation and administration of measurement data
- Logical design, self-explanatory
- Measurement protocol in accordance with standards
- Easy selection and configuration of common analytical scales
- English and German
- For automatic control of sieve shakers AS 200 control, AS 300 control, AS 450 control, AS 400 control, AS 200 tap, AS 200 jet
- For Windows® 2000 / XP / Vista / 7 / 8



Code	Description	Packaging
LLG09739250	Software EasySieve® for Sieve Shakers AS control / AS 200 tap / AS jet, Type Adapter RS232-USB	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

