

**Tripod stand****ISOLAB****Tripod stand**

Chrome plated steel. Tripods are meant to be used in conjunction with wire gauze, triangles or a sand bath. Holding ring with 150 mm diameter can be reduced to 125 mm with the supplied reduction attachment.



| Code        | Description                                                  | Packaging |
|-------------|--------------------------------------------------------------|-----------|
| LLG04008512 | Tripod stand, Diam. 150 mm, Adjustable height 200 ... 300 mm | 1 pz.     |
| LLG06242704 | Tripod stand, Diam. 125 mm, Total height 200 mm              | 1 pz.     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

