



## Ethidium Bromide Extractor

### Ethidium Bromide Extractor

For the rapid and safe removal of ethidium bromide from gel-staining solutions. Contains an activated carbon matrix, which quickly and easily removes >99 % of ethidium bromide. Each device can decontaminate up to 10 litres of gel-staining solution. After filtration, the decontaminated solution can be safely poured down the laboratory drain. The extractor funnel device fits most standard laboratory flasks and bottles (33 to 45 mm necks), and the unit includes a cap for storage between uses. Pack of 2



Code	Description	Packaging
LLG06282288	Ethidium Bromide Extractor, Type Ethidium Bromide Extractor	2 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

