

**Disinfectant for temperature-sensitive materials, Helipur H plus N****B. Braun Deutschland****Disinfectant for temperature-sensitive materials, Helipur® H plus N**

Liquid concentrate for disinfecting flexible endoscopes, anaesthetic accessories, rubber and plastic items.

- Also suitable for disinfecting metal, glass and ceramics
 - Effective against bacteria incl. T.B.-B, spores, fungi, lipophilic viruses including Hepatitis-B and HIV, coxsackie, adeno and vaccinia viruses
 - Can also be used in ultrasonic baths
 - Formaldehyde free
 - Material compatability reports are available from leading manufacturers.
- Use disinfectants safely. Always read the label and product information before use.

This product is a disinfectant containing biocides. Please inquire its availability and legal information.

Application concentration: 1.5% reaction time 15mins. = 15ml in 1L water

1.0 % reaction time 60 mins. = 10 ml in 1 L water



Code	Description	Packaging
LLG09191277	Disinfectant for temperature-sensitive materials, Helipur® H plus N, Type Dosing bottle, 1000 ml	1 pz.
LLG09191278	Disinfectant for temperature-sensitive materials, Helipur® H plus N, Type Jerrycan, 5000 ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

