

**Mixers MICROTRON MB 550/MB 800****KINEMATICA****Mixers MICROTRON® MB 550/MB 800**

The systems MICROTRON® MB 550 and MB 800 are specially designed for mixing, dispersing, homogenizing, emulsifying and size reduction of cuttable materials. Four different attachments made from glass with volumes of 125 ml, 250 ml, 500 ml and 1000 ml are available for the MICROTRON® MB 550. Attachments made from stainless steel with volumes of 2000 ml and 4000 ml are available for the MICROTRON® MB 800. Both systems are safety laboratory devices: grounded, double safety isolated, radio screened. A safety system prevents uncontrolled starting and therewith personal injuries: the mixer can only be started with well attached safety hood or lid. In case the hood or lid will be removed during operation the device will be stopped automatically. Both systems are designed to withstand high loads; they are space saving and always ready to use.

Mixing jars and protective hoods - please order separately. Mixer MICROTRON® MB 550 can only be started with attached safety hood.

Specifications**Motor input power****MB 550:****MB 800:****550W****800W****Speed range****MB 550:****MB 800:****Stepless, 1000 to 15,000rpm****Stepless, 1000 to 17,000rpm****Dimensions (Ø x H): 190 x 245mm**

Code	Description	Packaging
LLG06241879	Laboratory mixer MICROTRON® MB 550 230 V	1 pz
LLG06259416	Laboratory mixer MICROTRON® MB 800 230 V	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

