

**Cryogenic box Work2Store, PP****Heathrow Scientific****Cryogenic box Work2Storeä, PP**

Reduce sample handling time and the need for several racks when working with samples for cold temperature storage. Patent-pending opening and closing mechanism allows the rack to expand, creating needed room for sample manipulation. Then simply collapse the rack back to its compact form and it is ready for the freezer. Temperature range -80 to 121°C.

- Clear lid locks the rack together for safe handling and transportation
- Wells have a center depression to stabilize conical tubes and a cryogenic vial lock for one handed cap removal
- Stackable, autoclavable and with alphanumeric well identification
- Will fit into standard 2" (52mm) freezer racking



Code	Description	Packaging
LLG09194001	Cryogenic box Work2Store, PP, Type Work2Store, For tubes 1.5/2.0 ml	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

