

**Cryo-Safe Coolers****Bel-Art Products****Cryo-Safeä Coolers**

Use to protect samples on the bench, during transport, or in freezers with defrost cycle and fluctuating temperatures

Polycarbonate. A non-toxic freezing gel is isolated within the walls of the cooler, guaranteeing that labels or marks on the vials or tube are not damaged.

-The indexed and transparent (9.401 090, 9.401 091, 9.401 092) or gel-filled (9.401 093, 9.401 094) polycarbonate lid shows tube location and matches locations printed on the cooler

-Side handles (9.401 090, 9.401 091) or a wire handle for easy carrying

-Coolers can be stacked to save space and maximize freezer volume

-Rubber feet to prevent slipping

-All polycarbonate cooler body easy to clean

Store at -20°C for 24 hours prior to use for -15°C-Cryo-Safe-Boxes

Store at -20°C for a few hours until the temperature of the cooler falls below 0°C for -0°C-Cryo-Safe-Boxes



Code	Description	Packaging
LLG09401090	Cryo-Safe Coolers, For tubes 0.5 ml, 1.5 ml ml, dim. (l x w x h) 145 x 100 x 95 mm, Array 12	1 pz.
LLG09401091	Cryo-Safe Coolers, For tubes 0.5 ml, 1.5 ml, 2.0 ml ml, dim. (l x w x h) 145 x 100 x 95 mm, Array 12	1 pz.
LLG09401092	Cryo-Safe Coolers, For tubes 15 ml ml, dim. (l x w x h) 197 x 140 x 190 mm, Array 12	1 pz.
LLG09401093	Cryo-Safe Coolers, For tubes 0.5 ml, 1.5 ml ml, dim. (l x w x h) 243 x 157 x 146 mm, Array 32	1 pz.
LLG09401094	Cryo-Safe Coolers, For tubes 0.5 ml, 1.5 ml, 2.0 ml ml, dim. (l x w x h) 243 x 157 x 146 mm, Array 32	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

