

**Maxima H Minus First Strand cDNA Synthesis Kit with dsDNase****Thermo Scientific**

Sistema completo per la sintesi altamente efficiente del primo filamento di cDNA - fase di rimozione di DNA genomico integrata - temperatura di reazione incrementata: la sintesi può avvenire fra 42 e 65°C - alta resa con RNA fino a 20 kb - priming flessibile: oligo(dT)18, esameri random primer gene-specifici - applicazioni: sintesi del primo filamento di cDNA per RT-PCR, costruzione di librerie di cDNA, generazione di sonde per ibridazione, sintesi di RNA antisenso - include: mix enzimatico Maxima H Minus, Oligo(dT)18 e esameri random, RT Buffer 5X, dNTP Mix, acqua esente da nucleasi, dsDNase e dsDNase Buffer 10X.



Code	Description	Packaging
FL50K1681	Maxima H Minus First Strand cDNA Synthesis Kit with dsDNase	20 reaz.
FL50K1682	Maxima H Minus First Strand cDNA Synthesis Kit with dsDNase	100 reaz.



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

