

**FastDigest Aat II****Thermo Scientific**

5'...G A C G T#C...3' 3'...C#T G C A G...5' L'enzima di restrizione FastDigest AatII riconosce i siti GACGT<sup>A</sup>C e taglia alla temperatura di 37°C in 5-15 minuti, utilizzando il tampone universale FastDigest Buffer (isoschizomeri: Zral).



| Code      | Description        | Packaging  |
|-----------|--------------------|------------|
| FL5FD0994 | Fast Digest Aat II | 50 RXN (*) |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

