

**FastDigest Aat II****Thermo Scientific**

5'...G A C G T#C...3' 3'...C#T G C A G...5' L'enzima di restrizione FastDigest AatII riconosce i siti GACGT[^]C e taglia alla temperatura di 37°C in 5-15 minuti, utilizzando il tampone universale FastDigest Buffer (isoschizomeri: ZraI).



Code	Description	Packaging
FL5FD0994	Fast Digest Aat II	50 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

