



FastDigest Ajul

Thermo Scientific

5'...# 7(N) G A A (N)7 T T G G (N)11#...3' 3'...#12(N) C T T (N)7 A A C C (N)6 #...5' L'enzima di restrizione Ajul riconosce i siti $\wedge^{(7/12)}\text{GAA(N)7TTGG(11/6)}\wedge$ e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer.



Code	Description	Packaging
FL5FD1954	FastDigest Ajul	20 rxn (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

