

**FastDigest BamHI****Thermo Scientific**

5'...G#G A T C C...3' 3'...C C T A G#G...5' L'enzima di restrizione BamHI riconosce i siti G[^]GATCC e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer.



Code	Description	Packaging
FL5FD0054	FastDigest BamHI	800 RXN (*)
FL5FD0055	Fast Digest BamHI	2500 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

