



## FastDigest BglI

Thermo Scientific

5'...G C C N N N N#N G G C...3' 3'...C G G N#N N N N C C G...5' L'enzima di restrizione BglI riconosce i siti GCCNNNN<sup>^</sup>NGGC e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer.



| Code      | Description     | Packaging   |
|-----------|-----------------|-------------|
| FL5FD0074 | FastDigest BglI | 200 RXN (*) |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

