

**FastDigest BglII****Thermo Scientific**

5'...A#G A T C T...3' 3'...T C T A G#A...5' L'enzima di restrizione BglII riconosce i siti A[^]GATCT e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer.



Code	Description	Packaging
FL5FD0083	FastDigest BglII	100 RXN (*)
FL5FD0084	FastDigest BglII	200 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

