

**FastDigest BplI****Thermo Scientific**

5'...# 8(N) G A G (N)5 C T C (N)13#...3' 3'...#13(N) C T C (N)5 G A G (N)8 #...5' L'enzima di restrizione BplI riconosce i siti $^{(8/13)}\text{GAG(N)5CTC(13/8)}$ e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer.



Code	Description	Packaging
FL5FD1314	Fast Digest BpII	20 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

