

**FastDigest Bpu10I****Thermo Scientific**

5'...C C#T N A G C...3' 3'...G G A N T#C G...5' L'enzima di restrizione Bpu10I riconosce i siti CCTNAGC(-5/-2)^ e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer.



| Code      | Description        | Packaging  |
|-----------|--------------------|------------|
| FL5FD1184 | Fast Digest Bpu10I | 20 RXN (*) |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

