

**FastDigest BseGI****Thermo Scientific**

5'...G G A T G N N#...3' 3'...C C T A C#N N ...5' L'enzima di restrizione BseGI riconosce i siti GGATG(2/0)^ a 55°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: BstF5I, BtsCI, FokI).



Code	Description	Packaging
FL5FD0874	FastDigest BseGI	100 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

