

**FastDigest BseGI****Thermo Scientific**

5'...A C T G G N#...3' 3'...T G A C#C N ...5' L'enzima di restrizione BseNI riconosce i siti ACTGG(1/-1)^ e taglia a 65°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: Bse1I, Bsrl, BsrSI).



Code	Description	Packaging
FL5FD0884	FastDigest BseNI	100 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

