

**FastDigest BseXI****Thermo Scientific**

5'...G C A G C (N)8 #...3' 3'...C G T C G (N)12#...5' L'enzima di restrizione BseXI riconosce i siti GCAGC(8/12)^ e taglia a 65°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: BbvI, BstV1I).



Code	Description	Packaging
FL5FD1454	Fast Digest BseXI	20 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

