

**FastDigest Bsp120 I****Thermo Scientific**

5'...G#G G C C C...3' 3'...C C C G G#G...5' L'enzima di restrizione Bsp120I riconosce i siti G[^]GGCCC e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: PspOMI).



Code	Description	Packaging
FL5FD0134	FastDigest Bsp120 I	200 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

