

**FastDigest Bsp1407I****Thermo Scientific**

5'...T#G T A C A...3' 3'...A C A T G#T...5' L'enzima di restrizione Bsp1407I riconosce T[^]GTACA e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: BsrGI, BstAUI).



Code	Description	Packaging
FL5FD0933	FastDigest Bsp1407I	50 RXN (*)
FL5FD0934	FastDigest Bsp1407I	100 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

