

**FastDigest Eam 1105I****Thermo Scientific**

5'...G A C N N N#N N G T C...3' 3'...C T G N N#N N N C A G...5' L'enzima di restrizione Eam1105I riconosce i siti GACNNN[^]NNGTC e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: AhdI, AspEI, BmeRI, Dril, EclHKI).



Code	Description	Packaging
FL5FD0244	FastDigest Eam 1105I	100 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

