

**FastDigest EcoRI****Thermo Scientific**

5'...G#A A T T C...3' 3'...C T T A A#G...5' L'enzima di restrizione EcoRI riconosce i siti G[^]AATTC e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer.



Code	Description	Packaging
FL5FD0274	FastDigest EcoRI	800 RXN (*)
FL5FD0275	FastDigest EcoRI	2500 reaz.



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

