



FastDigest Fok I

Thermo Scientific

5'...G G A T G (N)⁹...3' 3'...C C T A C (N)¹³...5' L'enzima di restrizione FokI riconosce i siti GGATG(N)⁹ e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer.



Code	Description	Packaging
FL5FD2144	FastDigest Fok I	100 ul (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

