



FastDigest HhaI

Thermo Scientific

5'...G C G#C...3' 3'...C#G C G...5' L'enzima di restrizione HhaI riconosce i siti GCG^C e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: AspLEI, BstHII, CfoI, HinfI, HspAI).



Code	Description	Packaging
FL5FD1854	FastDigest Hhal	200 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

