

**FastDigest HindIII****Thermo Scientific**

5'...A#A G C T T...3' 3'...T T C G A#A...5' L'enzima di restrizione HindIII riconosce i siti A[^]AGCTT e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer.



Code	Description	Packaging
FL5FD0504	FastDigest HindIII	800 RXN (*)
FL5FD0505	Fast Digest HindIII	2500 reaz.



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

