

**FastDigest HpyF10VI****Thermo Scientific**

5'...G C N N N N N#N N G C...3' 3'...C G N N#N N N N N C G...5' L'enzima di restrizione HpyF10VI riconosce i siti GCNNNNN[^]NNGC e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: BstMWI, MwoI).



Code	Description	Packaging
FL5FD1734	FastDigest HpyF10VI	20 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

