

**FastDigest Mau BI****Thermo Scientific**

5'...C G#C G C G C G...3' 3'...G C G C G C#G C...5' L'enzima di restrizione MauBI riconosce i siti $CG^{\wedge}CGCGCG$ e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer.



Code	Description	Packaging
FL5FD2084	FastDigest Mau BI	20 UL (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

