

**FastDigest MboII****Thermo Scientific**

5'...G A A G A (N)8#...3' 3'...C T T C T (N)7#...5' L'enzima di restrizione MboII riconosce i siti GAAGA(8/7)^ e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer.



Code	Description	Packaging
FL5FD0824	FastDigest MboII	50 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

