

**FastDigest NcoI****Thermo Scientific**

5'...C#C A T G G...3' 3'...G G T A C#C...5' L'enzima di restrizione NcoI riconosce i siti C[^]CATGG e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: Bsp19I).



Code	Description	Packaging
FL5FD0573	FastDigest NcoI	20 RXN (*)
FL5FD0574	FastDigest NcoI	100 RXN (*)
FL5FD0575	Fast Digest NcoI	300 reaz.



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

