

**FastDigest NdeI****Thermo Scientific**

5'...C A#T A T G...3' 3'...G T A T#A C...5' L'enzima di restrizione NdeI riconosce i siti CA[^]TATG e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: FauNDI).



Code	Description	Packaging
FL5FD0583	FastDigest NdeI	100 RXN (*)
FL5FD0584	FastDigest NdeI	300 RXN (*)
FL5FD0585	FastDigest NdeI	1000 RXN. (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

