

**FastDigest Hin1II****Thermo Scientific**

5'... C A T G#...3' 3'...#G T A C ...5' L'enzima di restrizione Hin1II riconosce i siti CATG[^] e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: CviAII, FaeI, FatI, Hsp92II).



Code	Description	Packaging
FL5FD1834	FastDigest Hin1II	100 rxn (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

