



FastDigest NotI

Thermo Scientific

5'...G C#G G C C G C...3' 3'...C G C C G G#C G...5' L'enzima di restrizione NotI riconosce i siti GC[^]GGCCGC e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: CciNI).



Code	Description	Packaging
FL5FD0593	FastDigest NotI	20 RXN (*)
FL5FD0594	FastDigest NotI	50 RXN (*)
FL5FD0595	FastDigest NotI	150 reaz.
FL5FD0596	FastDigest NotI	250ul



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

