

**FastDigest RsaI****Thermo Scientific**

5'...G T#A C...3' 3'...C A#T G...5' L'enzima di restrizione RsaI riconosce i siti GT[^]AC e taglia a 37°C in 5-15 minuti usando il tampone universale FastDigest Buffer (Isoschizomeri: AfaI, CviQI, RsaNI).



Code	Description	Packaging
FL5FD1124	FastDigest RsaI	100 RXN (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

